



Email :  
sales@laboratoryinstrumentsmanufacturers.com

**Product Name :**  
Steep Plane With the Ruler and Accessories

**Product Code :**  
LIM-CAT-L0018-00012



**Description :**

Steep Plane With the Ruler and Accessories

**Technical Specification :**

Steep plane with the ruler and accessories

The structure is made of metal, so it provides resistivity and robustness.

It consists of 2 metal plates with a hinge at one end; one of the plates is used as a base, and the other one as an oblique plane.

The plane, that acts as a slope, represents graduated scale 20-40

cm (divided by 1 cm) on one side, so the distance could be measured along. It can be tilted and adjusted at any angle between 0 and 40°. Vertical scale 0-25 cm is also printed so the height of the slope could be measured and recorded. The instrument is completed with a plastic trolley with a hook, set of weights (including 1 x 100 g, 1 x 50 g and 1 x 20 g) and a container with a light string.

Steep Plane With the Ruler and Accessories, Steep Plane With the Ruler and Accessories Bulk Suppliers, Steep Plane With the Ruler and Accessories Tools, Steep Plane With the Ruler and Accessories Steep Plane With the Ruler and Accessories, Steep Plane With the Ruler and Accessories Manufacturers, Steep Plane With the Ruler and Accessories Suppliers from India, China, Kenya.

---

# Laboratory instruments manufacturers India

